RUTILIGHT

ES	Instrucciones de uso fuente de alimentación	3
EN	Instructions for use power supply unit	7
FR	Mode d'emploi Alimentation électrique	11
DE	Gebrauchsanweisung Netzteil	15
PT	Unidade de Alimentação power supply unit	19
IT	Istruzioni per l'uso Alimentatore elettrico	23

VALIDEZ

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RUTILIGHT

DESCRIPTION:

The RUTILIGHT® power supply is specifically designed to be used with the RUTILIGHT® device. Both allow lighting the surgical field in a comfortable and practical way, regardless of environment lighting.

On the front it has three separate outputs enabling to power up to 3 units of the RUTILIGHT® device simultaneously.

Each output has a switch with an indicator light. When an output is connected, the green indicator light on the switch turns on allowing to easily see that the RUTILIGHT® device is being powered. Otherwise, the green indicator light is off.

The outputs provide a DC voltage of 6.2V

The dimensions of the power supply are designed to easily fit anywhere but it features an extra weight and 4 rubber feet to guarantee good stability. The housing of the power supply is made of stainless steel.

The RUTILIGHT® power supply unit only uses RF energy for its internal function. Therefore, the RF emissions are very low and not likely to cause interference in nearby electronic equipment.


SYMBOLS USED

	Meets the requirements of Regulation 2017/745/EEC on medical devices
	Please read the operating instructions
SN	Serial No. Each power supply unit has a unique serial No. (20XX-XXXX). The first 4 digits indicate the year of manufacture.
	Class II equipment
	Identification of electrical and electronic equipment in accordance with Directive 2012/19/EC. <i>The product should not be disposed of with conventional household waste.</i>
	BF applied part

Switches (3)

Rutilight® connectors (3)

AC Mains connector

Mains ON/OFF switch

SPECIFICATIONS

The RUTILIGHT® power supply unit must be connected to an AC 100-240 V power source using a 4m long power cord without ground connection, which must be connected to the AC Mains connector located on the back panel of the power supply, next to the Mains ON/OFF switch.

The front panel of the power supply features the three connectors which allow to power up to 3 units of the RUTILIGHT® device simultaneously. Each connector has its own ON/OFF switch.

POWERING THE POWER SUPPLY UNIT

Turn on the power supply unit using the Mains ON/OFF switch located on the back panel.

Place the RUTILIGHT® device on the surgical instrument.

Connect the RUTILIGHT® device to the front panel connector(s) and when light is required in one of the devices, turn the corresponding switch to the up position. The (green) indicator light of the switch will come on.

To switch off the RUTILIGHT® device, turn the corresponding switch to the down position. The indicator light will go off.

It is recommended to turn off the Mains ON/OFF switch located on the back panel after using the RUTILIGHT® device(s).

Input voltage	100-240 V (50/60 Hz)
Output voltage	6,2 V DC +/- 5%
Output current	190 mA +/- 10% each output
Dimensions	44*200*200 mm
Weight	3 Kg
Accessories	- 4 m Power Cord (NON-GROUNDED POWER CORD EU IEC320 C7 H05VV) - RUTILIGHT® device – Part No. RL-01 Other power cords which can be provided by the manufacturer: - 3 m Power Cord (NON-GROUNDED POWER CORD USA IEC320 C7 NEMA 1-15P) - 3 m Power Cord (NON-GROUNDED POWER CORD AUS IEC320 C7 H05VVH2-F) - 3 m Power Cord (NON-GROUNDED POWER CORD UK IEC320 C7 H05VVH2-F) - 1,8 m Power Cord (NON-GROUNDED POWER CORD JPN IEC320 C7 JIS 8303)
Operating temperature	The RUTILIGHT® devices are connected to the power supply unit. The RUTILIGHT® devices are considered as an applied part of the assembly. However, they are not intended to be in contact with the patient except accidentally for a short time when the surgeon is operating. Therefore, we consider that the potential contact time with the patient of the applied part should not exceed one minute. Inside the operating room, the usual work environment is 18-25°C. In this environment, the device could reach a maximum temperature of 41-48°C. The device head could reach 56°C at maximum operating temperature

CLEANING RECOMMENDATIONS

The power supply unit must be disconnected from the mains AC source before cleaning the housing.

Only the outside of the power supply unit, a stainless steel housing, can be cleaned. The housing must be cleaned with a soft cloth dampened in a solution of neutral soap and water or with a disinfectant (such as quaternary ammonium or an alcohol disinfectant). The cleaning fluid must never be applied directly to the power supply unit.

However, the instructions of the cleaning agent manufacturer regarding material compatibility, dosage and application time shall be observed in all cases. Before using the power supply unit again, it must be completely dry.

SERVICING/REPAIR

Replacing or modifying parts on the unit is strictly forbidden and the power supply unit can only be used with the power cord provided and with a RUTILIGHT® device (part applied to the patient).

WARNING: The use of a power cord different from the one supplied could cause increased emissions or decreased immunity of the equipment. Repair or servicing of the equipment by the user is forbidden. In case of malfunction, please contact your authorized distributor. Only the repair service has access to the detail of all the parts (lists of components, descriptions and instructions) to ensure the correct repair of the power supply unit.

ENVIRONMENTAL CONDITIONS

STORAGE TEMPERATURE	-10°C to 60°C with 20% to 90% RH (non-condensing)
OPERATING TEMPERATURE	10°C to 30°C with 30% to 80% RH (non-condensing)
ATMOSPHERIC TEMPERATURE	620hPa to 1013 hPa
MAXIMUM OPERATING ALTITUDE	4000 m above sea level

The power source used by the RUTILIGHT® power supply unit should be that of a typical hospital environment. However, the RUTILIGHT® power supply unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

MANUFACTURER'S DECLARATION REGARDING ELECTROMAGNETIC COMPATIBILITY

The FA-RL01 power supply has been developed and tested based on the UN EN 60601-1-2 standard.

This product complies with the limit values established within the framework of this standard for the emission of and resistance to electromagnetic interference. Therefore, no interference with other equipment caused by the FA-RL01 power supply unit should be expected. However, it is possible that other units and/or accessories present in the operating environment do not comply with these limit values, causing interference with the FA-RL01 power supply unit, which could lead to malfunction.

Therefore, the following should be taken into account:

- all units and/or accessories in the operating environment should have a medical engineering approval;
- all equipment must be identified with the CE marking and have a valid declaration of conformity;
- units and accessories that do not meet the limit values or electromagnetic emissions must be duly identified.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONSEN 60601-1-2, CHAPTER 5.2.2.1, TABLE 1

RUTILIGHT® power supply unit is intended for use in an electromagnetic environment specified below. The customer or user of RUTILIGHT® power supply unit should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	RUTILIGHT® power supply unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RUTILIGHT® power supply unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / Flickers emissions IEC 61000-3-3	Fulfill	

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY EN 60601-1-2, CHAPTER 5.2.2.1, TABLE 2

RUTILIGHT® power supply unit is intended for use in an electromagnetic environment specified below. The customer or user of RUTILIGHT® power supply unit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV por contact ±2 kV, ±4 kV, ±8 kV, ±15 kV Air	±8 kV por contact ±2 kV, ±4 kV, ±8 kV, ±15 kV Air	Floors should be Wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ bursts IEC 61000-4-4	±2 kV for Power supply lines ±1 kV for I/O lines (input/output)	±2 kV for Power supply lines Not Applicable	Mains power quality should be that of a typical commercial and/or hospital environment.
Surges IEC 61000-4-5	±0,5 kV, ±1 kV; differential mode	±0,5 kV, ±1 kV differential mode	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% Ut (>95% dip in Ut) for 0,5 cycles 0% Ut (95% dip en Ut) for 1 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5s	<5% Ut (>95% dip in Ut) for 0,5 cycles 0% Ut (95% dip en Ut) for 1 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5s	Mains power quality should be that of a typical commercial and/or hospital environment. Mains power quality should be controlled, protected or free from these events. If continued operation during power mains interruption is required, RUTILIGHT® power supply unit should be plugged into a source of uninterruptible power supply or a battery.
Frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: Ut is the AC mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY **EN 60601-1-2, CHAPTER 5.2.2.1, TABLA 4**

RUTILIGHT® power supply unit is intended for use in an electromagnetic environment specified below. The customer or user of the RUTILIGHT® power supply unit should assure that it is used in such an environment

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
RF conducida IEC 61000-4-6 RF radiada IEC 61000-4-3	3 Vrms 150 kHz a 80 MHz 6 Vrms 150 kHz a 80 MHz 3 V/m 80 MHz a 2,7 GHz	3 Vrms 150 kHz a 80 MHz 6 Vrms 150 kHz a 80 MHz 3 V/m 80 MHz a 2,7 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the RUTILIGHT® power supply unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,17 \sqrt{P}$ 150 kHz to 80 MHz $d = 1,17\sqrt{P}$ 80 MHz to 800 MHz $d = 2,33\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^a . Interference may occur in the vicinity of the equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RUTILIGHT® power supply unit is used exceeds the applicable RF compliance level above, the RUTILIGHT® power supply unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the RUTILIGHT® power supply unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCE BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE RUTILIGHT® POWER SUPPLY UNIT
EN 60601-1-2, CHAPTER 5.2.2.1, TABLE 6

Rutilight® power supply unit is intended for use in a electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Rutilight® power supply unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the Rutilight® power supply unit as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 MHz to 80 MHz $d=1,17\sqrt{P}$	800 MHz to 800 MHz $d=1,17\sqrt{P}$	800 MHz to 1 GHz $d=2,33\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,74
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 to 80 MHz and 800 MHz, the separation distance is applied in the highest frequency range.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people

WARNINGS

The use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in faulty operation. If such use is necessary, this and any other equipment should be checked to ensure operational safety.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used maintaining a minimum distance of 30 cm (12 inches) from any part of the supply power unit (FA-RL01), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The power supply is disconnected by removing the plug. It is advisable not to place the power supply plug in such a way that it is difficult to disconnect it from the power source.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

ENVIRONMENTAL PROTECTION INFORMATION


Packaging

The packaging can be discarded with conventional household or paper/cardboard waste.

To protect the environment during use

In the case of a prolonged pause, it is advisable to switch off the unit to ensure operational safety and to achieve greater power saving.

Disposal

Identification of electrical and electronic equipment in accordance with Directive 2012/19/EC (waste from electrical and electronic equipment WEEE).
This symbol on the product or the packaging indicates that the product should not be disposed of with conventional household waste.

VALIDITY

These instructions for use supersede all previous versions.

TRADEMARK

RUTILIGHT® is a registered trademark of the Germans Trias i Pujol Foundation for Research on Health Sciences.

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