DESCRIPCIÓN

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RUTILIGHT


DESCRIPTION

The RUTILIGHT® device is an adaptable, versatile and portable lighting device with LED technology that can be attached to any surgical instrument, such as the handle of an electrosurgical knife.

It consists of a LED ring and a 3,2 m extensible cable.

The body of the RUTILIGHT device consists of medical-grade silicone and a medical grade PVC cable. The components of the RUTILIGHT® device (RL-01) are shown in figure 1.

Figure 1: RUTILIGHT®

The RUTILIGHT® device can only be connected to the FA-RL01 power supply (see figure 2).

RUTILIGHT® should be used in typical hospital environment.

The maximum operating voltage of the moving part of the device is 6.2V (+/- 5%).

ACCESSORIES

It is necessary to use the specific FA-RL01 power supply provided.

Figure 2: Front of the power supply provided with RUTILIGHT® device

Power switch (3) RUTILIGHT® connectors (3)

Figure 3: Back of the power supply provided with RUTILIGHT®

Output for connecting the power supply to the mains Power supply mains switch

INTENDED USE

The RUTILIGHT® is a single-use sterile light device for removable attachment to an electro-surgical pencil or other device intended to illuminate the site of observation/intervention.

The RUTILIGHT® device enables illumination of the surgical field in a convenient and practical way, allowing access to poorly-lit areas even with the use of headlights.

This improves the surgeon's field of vision in any condition, particularly in deep or difficult-to-access incisions. It is especially suitable for all types of open surgery, always in an operating room under normal lighting conditions.

The RUTILIGHT® device easily attaches to most surgical instruments, without interfering with their operation. It can also be mounted from one side on Yankauer cannulas, needle holders, sealers, coagulators, or any instrument from 5 to 8 mm in diameter.

The intended user is the surgical team inside the operating room, usually a surgeon, assistant surgeon, surgical nurse, anaesthesiologist, surgical technician or scrub nurse.

It is a sterile, single use device.

It provides 300 minutes (5 hours) of bright light.

The RUTILIGHT® device flashes 3 times as a warning signal to indicate that its useful life will expire in 30 minutes.

Important: For interventions where the use of the RUTILIGHT® devices may exceed 5 hours, it is recommended to have 2 RUTILIGHT® devices in the operating room in order to maintain the continuity of the additional light provided by the RUTILIGHT® device.

RUTILIGHT® has been tested and conforms to electrical safety and EMC standards according to EN 60601-1:2006 + Corr.:2010 + A1:2013 + AC:2014 + A12:2014; EN 60601-1-6:2010 + A1:2015 and UNE-EN 60601-1-2:2015

INDICATIONS

The RUTILIGHT® is indicated for opened surgical interventions and procedures requiring the use of electrosurgical pencil or other device and so that the area of/in a patient's body, where cutting and/or coagulation takes place can be illuminated, without any restriction for the age of the patient.

The intended clinical benefits of the RUTILIGHT® are indirect benefits to the patient because it doesn't interact with patient body but improves the surgeon visibility, particularly in case of narrow fields. RUTILIGHT® device

improves the surgeon visibility and visual perception of delicate structures, facilitates the surgery, reduces the operation time and minimize complication.

CONTRAINDICATIONS

It is not recommended for use in endocranial interventions, since there are no clinical trials available in direct contact with the central nervous system or in cardiac or central circulatory system interventions.

It is advisable not to leave the RUTILIGHT® device on when it is not being used, so as not to reduce its useful life, and to avoid looking directly at the LEDs when the RUTILIGHT® is on, as it may cause glare. In case of glare, it is recommended to look away immediately and wait about 5 seconds before resuming surgery.

RUTILIGHT® is not intended to use with X-rays.

WARNINGS AND PRECAUTIONS

It is recommended not to use the device in case the sterile pouch is damaged before use.

After use, dispose of the device in accordance with local regulations.

When the RUTILIGHT® device is not in use, it is recommended to switch it off before leaving it on the instrument table to avoid glare, and to preserve its useful life. Although the LED technology employed produces cold-light, direct and prolonged contact of the lenses on a sterile surgical drape may cause its weakening.

If it is necessary to exchange the device repeatedly during a surgical procedure, it is recommended, in order not to damage it and guarantee its correct functioning, to remove the instruments as gently as possible, without forcing the silicone base or pulling the cable.

During the surgery, bodily fluid or debris may stick to the lights lens and may affect the intensity of the emitted light. The user may use as sponge, wipe or gauze with sterile water to wipe the surface of the light lens.

Do not leave the RUTILIGHT® on when not in use, in order to preserve its useful life.

Handle the product with care. Do not drop or throw the device. Do not use excessive force when using RUTILIGHT® as this may damage the device. If damaged, do not continue using.

The RUTILIGHT® devices are connected to the power supply unit where are considered as an applied part of the assembly. RUTILIGHT® devices are not intended to be in contact with the patient except minimal contact for a short time when the surgeon is operating, not exceed one minute.

Inside the operating room, the usual work environment is 18-25°C. Under these conditions, the device could reach a maximum temperature of 41-48 °C. The device head could reach 56 °C at maximum operating temperature.

The device is a single use device, which is intended to be used exclusively for a single patient during a single procedure. Do not reuse, re-sterilize, reprocess and / or repack the device, because it may compromise the physical, chemical and biological properties of the device that are critical to the performance and safety of the device. This may lead to a serious deterioration of a patient's, user's or other person's state of health due to e.g. recurrence and infection (including cross infection).

Consult the instruction manual of the power supply (FA-RL01) for use and EMC warnings.

The cable of RUTILIGHT® device contains Di(2-ethylhexyl) phthalate (DEHP). Based on animal studies, significant exposure to DEHP could interfere with the normal development of the male reproductive tract. When using RUTILIGHT® devices, the risks for children, pregnant and nursing women is evaluated as very low because of limited potential contact of the cable with patient fluid pathway.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

INSTRUCTIONS

1. Open the pouch and remove the RUTILIGHT® device with sterile gloves to place it on the surgical field.
2. Handle the product with care without excessive force.
3. Mount the device by inserting the surgical instrument that is required at that moment through the central hole until it is fixed and does not move.
4. When mounting the device on an electrosurgical knife, attach the RUTILIGHT® device to the end of the base, without overlapping on the active terminal (blade, needle or clamp). Ensure that it is perfectly in place in the handle with the cable on the opposite side of the push buttons.
5. Extend the cable and plug it into the power supply.
6. Press the main switch on the back of the power supply and, when the device is needed, turn on the front switch to which it is connected to start it.
7. After 270 minutes (4h30m) of use, the light flashes briefly to indicate that it will turn off in 30 minutes. If more lighting time is needed when the useful life of the RUTILIGHT® device has expired, it will be necessary to use a new device.

8. After using RUTILIGHT® device, it should be disposed of according to local regulations for disposable surgical instruments.

PACKAGING, STERILITY AND STORAGE


RUTILIGHT® is supplied in sterile individual units. The intact sterile pouch protects and maintains the device sterilized and ensures sterility until the expiry date, provided that the storage has been correct, in its original packaging in a cool, dry place and protected from direct sunlight.

The RUTILIGHT® device should not be used after the expiry date (see label).

Before using the RUTILIGHT® device, make sure that the packaging is not damaged. If it is, it may no longer be sterile and should not be used in any case.

The RUTILIGHT® device has been sterilized with ethylene oxide.

SYMBOLS USED ON THE PACKAGING

	Part No. (P/N)
	Sterilized by ethylene oxide
	Use before the expiry date
	Keep away from direct sunlight
	Temperature limits
	Not re-sterilize
	Read the instructions for use
	Complies with the requirements of Regulation 2017/745 on medical devices
	Manufacture date
	Lot number
	Keep dry
	Do not reuse
	Do not use if packaging is damaged
	Manufacturer

	Medical Device
	phthalates
	Device restricted to use by a physician
	Warning
	Unique Device identifier
	Applicable part

Each box contains 10 or 25 RUTILIGHT® device units (P/N RL-01), each one packaged in a sterile, disposable pouch.

VALIDITY

These instructions for use supersede all previous versions.

TRADEMARK

RUTILIGHT® is a registered trademark of the Germans Trias i Pujol Foundation for Research on Health Sciences.

DECLARATION OF CONFORMITY

The RUTILIGHT® device complies with the requirements of Regulation 2017/745 on medical devices. It bears the symbol CE 2797.

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